

HUBBARD COMMUNICATIONS OFFICE
Saint Hill Manor, East Grinstead, Sussex
HCO BULLETIN OF 24 AUGUST 1971
Issue II

Remimeo

ASSISTS ADDITION

(Refers to HCO B 23 July 71
which revised HCO B 15 Mar 71)

In running Version A and Version B of a Hold it Still temperature assist DO NOT RUN IT OVER OUT RUDS.

Often a pc is ill because his ruds in life are out.

The Temperature Assist will only fail if the pc had an In Life ARC Break, Present Time Problem or Withhold. Sometimes these alone will change the temperature for the better.

Each rud of course must be taken Earlier Similar to F/N if no F/N on the first answer given to a rud.

So Add to Version A C/S at the start:

0. Fly all ruds.

And to Version B C/S at the start:

0. Fly all ruds.

These are done with the pc holding the cans. The walkabout version is done off the cans but checked.

Overrun is rehabbed. If no rehab then flatten by further running.

When the pc off the cans in Version B has a cognition he should be put back on the cans and checked. Usually an F/N will be found.

It has not been fully determined how many times a pc can be run on "Hold it Still".

But if the pc has been run before on the process and does not run, it should be checked for overrun and rehabbed.

L. RON HUBBARD
Founder

LRH:sb.rd